

# Germline Variant Calling in Formalin-fixed Paraffin-embedded Tumours

Eva Yap, MSc. Student

October 31, 2017

# Overview

- 1 Background
- 2 Research Question
- 3 Methods
- 4 Results
- 5 Conclusions

# Overview

- 1 Background
- 2 Research Question
- 3 Methods
- 4 Results
- 5 Conclusions

# Germline variants have important clinical implications

## Cancer Predisposition


- Preventive measures
- Sibling testing




## Pharmacogenomics

- Treatment tolerance
- Adverse drug reaction (ADR)

### Patients with same diagnosis



# The tumour genome contains germline information



# Clinical tumour sequencing could be a practical, cost-effective approach to provide germline testing

**Patients undergoing clinical tumour sequencing**




**Patients with potential germline variants**



**Downstream confirmatory testing**



# Challenge: Distinguishing between germline and somatic variants in the tumour genome



# Overview

- 1 Background
- 2 Research Question**
- 3 Methods
- 4 Results
- 5 Conclusions



# Research question


Can we accurately identify germline variants in tumour genomes?

# Overview

- 1 Background
- 2 Research Question
- 3 Methods**
- 4 Results
- 5 Conclusions

# Variant allele frequency (VAF)

**Tumour content**



100%

**Sequencing**



Heterozygous variant  
VAF = 50%



Homozygous variant  
VAF = 100%


# VAF in tumour specimens can deviate from diploid zygosity

DNA damage induced by formalin (e.g. fragmentation and sequence artifacts)




# Somatic VAF in tumour specimens can deviate from diploid zygosity

Mixture of tumour and normal cells




# Somatic VAF in tumour specimens can deviate from diploid zygosity

Tumour heterogeneity



# Study Design



		Predicted variant status	
		Germline	Somatic
Detection in matched blood	Present	True positive	False negative
	Absent	False positive	True negative

$$\text{Sensitivity} = TP / (TP + FN)$$



$$\text{PPV} = TP / (TP + FP)$$

# Overview

- 1 Background
- 2 Research Question
- 3 Methods
- 4 Results**
- 5 Conclusions





# VAF of germline variants in blood and FFPE tumours



# Sensitivity in detection of germline variants at different VAF thresholds

VAF (%)	False Negative	True Positive	Sensitivity	95% CI
10	0	1981	1.0	1.0–1.0
15	13	1968	0.99	0.99–1.0
20	46	1935	0.98	0.97–0.98
25	77	1904	0.96	0.95–0.97
30	117	1864	0.94	0.93–0.95
35	192	1789	0.90	0.89–0.92
40	313	1668	0.84	0.83–0.86
45	458	1523	0.77	0.75–0.79

# VAF of germline and somatic variants in FFPE tumour



High positive predictive value can be achieved for referral of germline variants to downstream confirmatory testing

VAF (%)	False Positive	True Positive	Total Calls	PPV	95% CI
10	431	1981	2412	0.82	0.81–0.84
15	319	1968	2287	0.86	0.85–0.87
20	273	1935	2208	0.88	0.86–0.89
25	245	1904	2149	0.89	0.87–0.90
30	203	1864	2067	0.90	0.89–0.91
35	178	1789	1967	0.91	0.90–0.92
40	146	1668	1814	0.92	0.91–0.93
45	118	1523	1641	0.93	0.91–0.94

# Overview

- 1 Background
- 2 Research Question
- 3 Methods
- 4 Results
- 5 Conclusions**

# Conclusions

- ① A VAF approach demonstrates high sensitivity and precision at separating between germline and somatic variants.

# Conclusions

- ① A VAF approach demonstrates high sensitivity and precision at separating between germline and somatic variants.
- ② At 30% VAF threshold, sensitivity for detection of germline variants in FFPE tumour is 0.94 and the positive predictive value for referral to downstream confirmatory testing is 0.90.

# Conclusions

- ① A VAF approach demonstrates high sensitivity and precision at separating between germline and somatic variants.
- ② At 30% VAF threshold, sensitivity for detection of germline variants in FFPE tumour is 0.94 and the positive predictive value for referral to downstream confirmatory testing is 0.90.
- ③ Germline variants could be accurately identified in FFPE tumour sequencing.



# Acknowledgements

Dr. Aly Karsan  
Dr. Jennifer Grants  
Dr. Jeremy Parker  
Dr. Kieran O'Neill  
Dr. Marion van den Bosch  
Dr. Rawa Ibrahim  
Dr. Sergio Martinez-Hoyer  
Angela Mo  
Patrick Coulombe  
Rod Docking  
Jenny Li  
Deborah Deng  
Helen Lin

Committee members:

Dr. Martin Hirst  
Dr. Ryan Morin

Centre for Clinical Genomics

Canada's Michael Smith  
Genome Sciences Centre

Patients from The OncoPanel  
Pilot study (H14-01212)